

# PRODUCT SPECIFICATION

REF 726000 CEDERROTH EYE & WOUND CLEANSING SPRAY



## 1 GENERAL PRODUCT DESCRIPTION

### Description

A portable spray bottle filled with 150 mL physiological (0.9%) unbuffered sodium chloride solution, preservative free and sterile. The nozzle is suitable for use on both eyes and wounds. The actuator is protected by a plastic cap and a tamper seal. The nitrogen propellant is not in contact with the solution. The solution remains sterile until last drop. The Cederroth Eye & Wound Cleansing spray is intended to be used for rinsing and cleansing out dust and dirt from eyes and wounds.

### Storage information

Store in temperatures between 0°C and +40°C.



## 2 DETAILED PRODUCT INFORMATION

### Content per consumer unit

#### Material

Liquid	Sterile 0.9% Sodium Chloride (NaCl). Preservative free.
Can	Aluminum.
Bag-on-valve	4 layer PP-OPA-ALU-PET
Actuator	White POM
Cover CAP	Transparent PP

**Shelf life:** 5 years from production date

## 3 LOGISTIC INFORMATION

	CU	SKU	Pallet
<b>GTIN:</b>	7310617260005	7310807260006	7312027260002
<b>Content:</b>	150 ml	12 CU	24 SKU/layer / 5 layers
<b>Size (w x l x h mm):</b>	53 x 53 x 162	168 x 221 x 191	Pallet 800 x 1200x 1105
<b>Gross weight:</b>	191 g	2,5 kg	317 kg
<b>Packing material:</b>	-	Carton 140 g	-

